



# Sample QC: Absorbance or Fluorescence?

## Absorbance

Absorbance methods directly measure the amount of light absorbed by a sample at a specific wavelength. The absorbance of light is proportional to concentration of the sample, allowing quantification.

- |  |   |
|--|---|
| <p> <b>Microvolume</b></p> <ul style="list-style-type: none"> <li>• 1 <math>\mu</math>L sample volume</li> <li>• No sample dilutions required</li> </ul> | <p> <b>Cuvette</b></p> <ul style="list-style-type: none"> <li>• Larger sample volume</li> <li>• Higher sensitivity than microvolume</li> </ul> |
|--|---|

## Fluorescence


Fluorophores are bound to the sample of interest and excited with a specific wavelength of light. Emission is measured at a higher wavelength. Sample fluorescence is compared to known standards.

- 0.5 mL PCR Tubes**
- 200  $\mu$ L assay volume
  - 1 - 20  $\mu$ L of sample
  - Choose from a wide variety of fluorophores and assay kits




## DeNovix DS-Series Applications

### Absorbance Apps



### Fluorescence Apps



## Feature Comparison: DS-11 FX+

<p><b>0.75</b> ng/<math>\mu</math>L</p> <p>Not Analyte Specific</p> <p><b>37,500</b> ng/<math>\mu</math>L</p> <p><b>0.5 - 1 <math>\mu</math>L</b></p> <p>Detects Co-Extracted Contaminants</p> <p><b>Zero</b> Cost / Sample</p> <p>Load and Measure in Seconds</p>	<p><b>Sensitivity</b></p> <p><b>Specificity</b></p> <p><b>Dynamic Range</b></p> <p><b>Sample Volume</b></p> <p><b>Contamination Detection</b></p> <p><b>Assay Cost</b></p> <p><b>Speed</b></p>	<p><b>0.0005</b> ng/<math>\mu</math>L</p> <p>Highly Specific</p> <p><b>4,000</b> ng/<math>\mu</math>L</p> <p><b>1 - 20 <math>\mu</math>L</b></p> <p>No Contamination Information</p> <p>Range of Assays</p> <p>Assay Setup Required</p>
--	--	---

## Summary

### Absorbance

- Quick and easy method to quantify a wide range of samples without any assay cost or set up time.
- Possible over a wide dynamic range, covering all commonly used samples concentrations.
- Contaminants can be detected, allowing additional quality control.

### Fluorescence

- Enhanced sensitivity allows quantification to sub-picogram per micro liter concentrations.
- Assays are specific for the analyte under investigation, leading to highly specific quantification without interference from contaminants.
- Requires assay kits and known-concentration standards.

## Combining Measurements

Both UV-Vis Absorbance and Fluorescence have advantages depending on the application and sample being quantified. The **DeNovix DS-Series** offers several instrument models with combinations of microvolume absorbance, cuvette absorbance and fluorescence measurement modes.




**DS-11 FX+**  
Spectrophotometer / Fluorometer

1  $\mu$ L UV-Vis, Cuvette, and Fluorescence in One Instrument



**DS-8X+**  
Eight Channel Spectrophotometer

1  $\mu$ L UV-Vis (x8), Cuvette, and Optional Fluorometer Module



### Free Fluoro Assay Evaluation Kit

- ✓ dsDNA or RNA assays
- ✓ 50 samples
- ✓ Free Shipping to Your Lab

