

The background of the entire page is a dark blue gradient with dynamic, glowing light trails in various shades of blue. These trails are curved and flow across the frame, creating a sense of motion and energy. The light trails are most prominent in the upper right and lower right areas, while the left side is darker.

DeNovix<sup>®</sup>

Fluorescence Measurement  
**Best Practices**

## Introduction

Fluorometric quantitation assays for nucleic acids and protein are now routine in many life science research labs. As with many fluorescence techniques, it is important to apply best practices to ensure reliable results and reduce user error.

This document provides some best practice tips for making fluorescent measurements using the DeNovix DS-11 Series.



## LED Selection

Ensure that the excitation source (LED) and emission filters are appropriate for the fluorophore or assay of interest (Table 1). DeNovix Fluorescence instruments are equipped with four user-selectable LED light sources. See Figure 1 for a list of common fluorophores.

LEDs	Excitation Filter Range	Emission Filter Range
UV (375 nm)	361 – 389 nm	435 – 514 nm
BLUE (470 nm)	442 – 497 nm	514 – 567 nm
GREEN (525 nm)	490 – 558 nm	565 – 650 nm
RED (635 nm)	613 – 662 nm	664 – 740 nm

## Sample Tubes

- Use only thin-walled, clear 0.5 mL PCR tubes for sample measurements.
- Use 200  $\mu$ L final volumes for all standards and samples.
- Do not label the side of an assay tube, as this could interfere with the sample measurement.
- Ensure that the tube is clean and dry on the outside before inserting into the DS-11 Series tube chamber. Moisture and condensation on the tube surface may lead to measurement errors.


## Sample Preparation

- Protect dye reagents and working solutions from light.
- Confirm that fluorophore assay reagent is compatible with sample buffer components.
- Follow assay manufacturer's recommendations regarding standard curve concentrations.
- Treat all samples and standards identically in terms of volumes, incubation times and temperature.
- Measure all standards and samples within the assay manufacturer's recommended assay time frames.

## Sample Measurements

- Ensure that the assay tubes (and solutions) are at room temperature at the time of the measurement. Temperature fluctuations may impact the accuracy of the assay.
- Ensure that sample solutions are homogeneous and well-mixed before sampling. Avoid introducing air bubbles into the sample solution when mixing sample.
- Close the fluorometer chamber lid before tapping the Measure button.

## Figure 1: Common Fluorophores



## Fluorescence Quantification Video

Learn how and when Fluorescence Quantification can be the better method to quantify biomolecules when using the DS-11 Series or QFX Fluorometer.



Note: Fluorescence results are reported in relative fluorescence units (RFU). Differences in RFU between different fluorometers are expected. Use a standard curve for accurate quantitation.

*Qubit™, Quant-iT™, Alexa Fluor™, PicoGreen™, Oligreen™, RiboGreen™ and SYBR™ are registered trademarks of Thermo Fisher Scientific and its subsidiaries and are used for identification and references purposes only. DeNovix, DeNovix products and this document are not endorsed or authorized by or in any way affiliated with Thermo Fisher Scientific. QuantiFluor™ is a registered trademark of Promega Corporation and is used for identification and references purposes only. DeNovix, DeNovix products and this document are not endorsed or authorized by or in any way affiliated with Promega Corporation.*



## Contact

---

+1.302.442.6911  
[info@denovix.com](mailto:info@denovix.com)

DeNovix Inc.  
3411 Silverside Road  
Hanby Building, Suite #101  
Wilmington, DE 19810



## Connect

---

@denovix on all platforms

[denovix.com/community](https://denovix.com/community)

Use #denovixlove for a chance to be featured!



## Explore

---

[Absorbance and Fluorescence Quantification eBook](#)  
[TN 156: Absorbance and Fluorescence Quantification](#)  
[TN 162: dsDNA Assay Concentration Limits](#)  
[TN 200: DeNovix RNA Assay Performance Data](#)  
[Absorbance or Fluorescence Infographic](#)  
[DS-11 Series Spectrophotometer / Fluorometer Info](#)  
[QFX Fluorometer Info](#)